

Water Baths





WATER BATHS

Laboratory water baths provide accurate and stable temperature control for sample heating applications. Features include anti-dry protection, seamless stainless steel chambers, customizable openings, and efficient heat distribution for reliable laboratory performance.



ANTRI DRY TECHNOLOGY

SEAMLESS RESERVOIR

SPILL SAFE CONTROLLER

ADJUSTABLE RING LID

TECHNICAL SPECIFICATIONS

Model	WB100-1F	WB100-2F	WB100-4F	WB100-6F	WB100-8F
Temperature Range	RT - 99.9°C				
Capacity	3L	6L	12L	18L	24L
Power	300W	600W	800W	1200W	1500W
Accuracy	±0.1°C				
Working Chamber Size (mm)	160 x 160 x 130	300 x 150 x 140	300 x 300 x 140	450 x 300 x 130	620 x 300 x 120
Anti-Dry	Yes				



WATER BATH WBS SERIES

WBS Series water baths provide accurate and stable temperature control for laboratory heating applications. Features include stainless steel chambers, anti-dry protection, transparent covers, and uniform heat distribution. The WBS-6 Pro model includes a built-in magnetic stirrer for enhanced temperature consistency.



ANTI DRY TECHNOLOGY

TALL TRANSPARENT COVER

STAINLESS STEEL TANK

STIRRER OPTION

TECHNICAL SPECIFICATIONS

Model	WB-2S		WB-6S	WBS-6Pro
Speed / Increment	-			100-1200RPM / 50RPM
Capacity	2L	6L		6L
Temperature Range	RT-100°C			
Bath Size	Φ150 x 150 x 105mm(depth)	200 x 200 x 130mm(depth)		
		2L flask or 3L beaker		
Product Size (mm)	195 x 190 x 260	260 x 265 x 325		260 x 265 x 325
Display Accuracy	0.1°C			
Power	300W	800W		
Voltage	AC110/220V			
Timing Range	1min-23h59min			



OIL BATH WITH MAGNETIC STIRRER

Oil baths with magnetic stirrers provide accurate heating and mixing for laboratory applications. Suitable for water or oil baths, features include precise temperature control, over-temperature protection, stainless steel chambers, and efficient heat distribution.



DUAL USE: WATER & OIL

HIGH TEMP UP TO 250°C

STAINLESS STEEL BODY

INCL STAND & PROBE

TECHNICAL SPECIFICATIONS

Model	OB-2C	OB-2Pro	OB-5Pro	OB-8Pro
Control mode	Digital Display & Knob	Digital Display & Button & Control Knob		
Speed range	50-1500rpm		50-1200rpm	
Timing range	-	1min~12h	1min~12h	1min~12h
Setting temp	RT-250℃			
Max mixing volume	1L flask/2L beaker		2L flask/5L beaker	5L flask/10L beaker
Inner dia./Depth	150/80mm		200/105mm	250/135mm
Toggle forward/reverse	-	Switchable		
Temp probe	External PT1000			
Power	300W		500W	800W
Voltage	AC110/220V			
Dimension	285 x 180 x 180mm	285 x 180 x 180mm	355 x 240 x 215mm	415 x 300 x 245mm